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Version 3.2.1.21751 of DWOS is a maintenance release for the 3Series, 5Series and 7Series as well as standalone software. The update features the following improvements:

Scan related issues:

- Fixed bug causing the scan workflow of an XOrder with implants to loop.
- Fixed bug during the scanning procedure of a simple coping with waxing.
- Enabled occlusion key re-scan when scanning a model with positive antagonist.

CAD related issues:

- Fixed inverted surfaces of connector when using Adapt Anatomy function on a full pontic.
- Fixed bug causing an edited connector to be recomputed after validating shape modifications on a connected tooth.
- Improved computing of Mirror Anatomy function to avoid generating 2D surfaces when the extracted area includes an artifact.

Implant related issues:

- Computing of screw hole through overpressed restorations on custom abutments.
- Improved merging result of custom abutment with waxing.
- Fixed crashing issues when editing an order containing a bridge over custom abutments.
- Fixed bug causing certain implants interface to be computed incorrectly after merging surfaces of a custom abutment.

Also:

- Fixed issues causing imports of XOrders to fail.
 - Enabled viewing of attached files in Administrative Application.
 - Now integrated Case Manager connection within the DWOS Inbox.
-

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