Release Notes DWOS Lava Edition 3.9.1

Version 1.0 8/29/2014

Dental Wings team has made available a new upgrade of its DWOS Lava Edition software. Version 3.9.1.29487 is compatible with the 3M scan ST scanners. The following is the complete list of new features and improvements:

New features :

- Crown & Bridge design
 - Shaping Tools
 - Transforms on anatomies allow smoother deformation
 - Use of Automatic Adjust Contact Points without defining mesial and distal points
 - Egg shell option for diagnostic crowns
 - Shoulderless Telescopic Crown
 - Insertion axis of attachments can now be edited in groups and is linked to telescopic axes
- Custom Abutments
 - New prosthesis type “Custom Gingiva Former” (custom implant healing abutment/cap)
- Model Builder
 - Saw Cut models (compatible with Baumann and Dentona plates)
 - Zimo Models with DWOS occlusor
 - Mono block models
 - Improved tagging function : embossed, engraved, plate
 - Add stabilization bars
- Implant bars and bar pillar design
 - Bar profile Editor
 - New default profiles : Cross, Lambda, Trapezoidal and Staircase retention
 - Create, save and export your custom profiles
 - New Bar Design functionalities : chamfer plane, attachments and spacers blocks
 - Add cutouts on bars to allow for retention area while lightening the bar
- Wizard
 - Custom Abutment wizard
 - Partial Framework Design Wizard

Improvements:

- Crown & Bridge design
 - Shaping tools
 - Global/Local deformation mode on Clinical handles editing
 - Automatic synchronization of anatomy and overlay cervical lines
 - Adaptation environment improvements
 - Retention beads on all restorations
 - New tool compensation parameters for substructure restorations
- Custom Abutments
 - Virtual waxing
 - Edit connectors, parallelism
 - Gingiva base line edition
 - Add dental bands, retention beads and attachments
 - Wax-up
 - Add dental band, retention beads and attachments
 - Improvements on design workflow and performances
- Model Builder
 - New Linear workflow
- Implant bars and bar pillar design
 - Attachment hole compensation upon merging to respect minimum thickness
 - Improved initial proposal (follow the curve of initial arch)
- Wizard
 - Ability to take screenshot
 - Attach files to orders
 - Desktop icon
 - Add adjacent scanning area step in Partial in wizard mode
- Partial
 - Inferior bar : new expansion parameter
 - Smoother damping of clasps
- Others
 - Allow dual arch cases through Scan import
 - New Hyperdent open output
 - New CADesthetics output format

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