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Version 1.2 5/13/2015

Dental Wings team has made available a maintenance release of its DWOS software, version 5.0.1.30084. This upgrade is compatible with all scanners as well as standalone software licenses. The following is the complete list of improvements and bug fixes:

Improvements

- Introduction of a new bar profile : Munich bar (inverted T-shape)
- Implant bar : shoulder support on bar pillars
- New retention elements for hybrid and wrap-around bars
- Orthodontic archiving: possibility to define the height of the model
- Improved rendering and file size of 3D PDF reports
- Updated labels for Czech and Russian languages

Bug fixes

- Stabilization work on waxing merge process (waxing copy, waxing on stumps and waxing on implants)
- Fixed issues when routing cases with dental bands and prehension items in the multi-die station
- Fixed cervical handles behavior when adjusting contact points on pontics
- Fixed some import issues for orders with pontics or invalid/not present implant kits
- Stabilization work on custom abutments for several implant kits (interface connection and screw channel problems).
- Stabilization work to prevent multiple transaction errors when using attachments or defining margin line
- Fixed some issues with 3/4 crowns over abutments
- Fixed issues while exporting to Zeno output format 3.2 and 4.0

- Fixed cases in order creation where the implant kit selection was not registered
- Fixed merging problems in the partial scan workflow on a 5 Series
- Fixed an error appearing when trying to scan with the vertex holder
- Full Denture: fixed computation and merging problems in the full denture module regarding the gingiva

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