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Version 1.2 9/28/2016

Dental Wings team has made available a maintenance update of its DWOS software. Version 6.1.1 is compatible with all scanners as well as standalone software licenses. The following is the complete list of bug fixes:

Bug Fixes

Pre-position anatomy

- Fixed erratic computing of hygienic and saddle pontics
- Improved reliability of Clinical Handles editing
- General stabilization of the feature

Global Axis usability

- Outer insertion axis of custom abutments and telescopic copings being grouped by default, it is now possible to use the *Axis Group* function to customize the groups.
- There is additionally a new option in the insertion axis window in the scan station to do so.

Partial Framework design

- Fixed misleading display of undercut filling on blocked models to display actual blocking wax
- Fixed issues causing the computation of occlusal rest to fail
- Mesh computation optimization

Others

- Fixed issues that prevented to adjust contact points on orders incoming through *Scan import station*
- Fix a bug to be able to add more than one support in the Rapid Prototyping Manufacturing (RPM) station
- Fixed bug in the display of *Margin Line Verification* when setting the insertion axis in cases with multiple bridge orders
- Fixed issues when importing XOrders containing a 3D files in PLY format (cases coming from 3M True Definition Scanner)
- Improved reliability when sending cases to DWOS Connect (avoid network error)
- CAD Inbox maintenance

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