# Release Notes - Notes de mise à jour

## DWOS 6.4.0

Version 1.0 27/01/2017

---

Dental Wings team has made available a maintenance update of its DWOS software. Version 6.4.0 is compatible with all scanners as well as standalone software licenses. The following is the complete list of improvements and bug fixes:

### Enhancements

#### Bite Splint Design

- Integrated minimum material thickness management

#### Virtual Model Builder

- Possibility to create flat-bottomed removable gingiva
- Widen parameter boundaries to increase flexibility
- Completed rework of parameter visual description for model and stump types

#### Partial Framework

- Possibility to add positive and negative attachments to the partial framework designs to create for example, retention elements and holes

### Bug Fixes

#### Virtual Model Builder

- Fixed problems while generating virtual implant model designed with a removable gingiva
- Improve computation speed of virtual model designed with a removable gingiva
- Correct the antagonist positioning in the virtual model builder station

#### Crowns and Bridges

- Improve default connector proposal on wide but small posterior restorations
- Fixed exceptions when using Anatomy Kits Editor application
- Fixed a case when adjusting contact points would create artifacts inside the intrados.
- Inverted transparency order to allow seeing the overlay through the wax-up in virtual waxing tools designer
- Fixed a margin projection bug when creating crown over post & core

- Improved natural connector final mesh quality
- Improved designer stability of 360 dental bands
- Preventing the telescopic crown top spline to create intersections
- Improved tooth chain initial proposal for cases with very low vertical space

#### Implants

- Fixed a problem of closed angulated screw channel
- Fixed a problem when using maximum body of some compatible implant kits which were failing the design parameters validation
- Fixed hyperdent (hdpp) output format for implants and implant bridges
- Improved support for ti-base abutments when enabling cylindrical base option
- Fixed 3/4 crowns on implant overlay edition

#### Partial Framework

- Fixed a problem while adding material once the partial framework had been merged
- Survey axis now keeps its current position when reopening the survey designer
- Wax-up scan can now be hidden

#### Full Denture

- Fixed base plate visibility after merge
- Improved merge stability

#### DWOS Connect

- Shipping status is now properly synchronized and being sent back from the laboratory to the dentist

#### Inbox

- Several minor bug fixes to enhance the use of the inbox application view
- Re-Assign controller is now visible on virtual model builder cases
- Fixed import of subcontracted Model Builder orders
- Fixed crashes with the item preview panel
- Improved memory management

#### Others

- Set correct date when generating output zip file
- Unify wax-up nomenclature throughout the software
- Fixed numbering inconsistencies between pictures and parameters in dialogs

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