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Version 1.0 9/27/2017

---

Dental Wings team has made available a new upgrade of its DWOS software. Version 7.1.0 is compatible with the following Dental Wings scanners: 7Series, 3Series, iSeries as well as standalone software licenses. This version is not compatible with 5Series scanners.

The following is the complete list of new features and improvements:

### New features

#### Crown and Bridge / Implant

- Tooth Chain improvements: activate independently interproximal and occlusal contacts (deactivate occlusal when allowed space is too low for good adjustment)
- Mirror anatomy on natural tooth selection on diagnostic wax-up surface
- 3-point repositioning available on Straumann implant kits in Scan station
- Possibility to create bridges on second layer

#### Model Builder

- Model Builder model and stump parameters are now material driven. Each parameter can be customized per material
- Improvements of removable flat bottom gingiva
- Create analog socket protection when creating virtual model with implants
- Add attachment on virtual models: create gingiva retentions, control window, lingual base

#### Partial Framework

- New tool available called “spacers” to simulate wax spacing using preparation wax in traditional method
- Add retention beads on partial framework

### Improvements

#### Crown and Bridge / Implant

- Undo/redo support while editing connectors
- Side view when editing telescopic crown
- Show/hide implant interface while drawing margin line

## Model Builder

- Zimo and Witlaken occlusor : customizable pin diameter for fine tuning while printing models

## Partial Framework

- High resolution mode for waxing tools available in design mode and merged state
- Undercut designer : improved axis definition view with rotation angle information displayed

## Others

- Undo/redo support while adding attachments
- Improved translations for French, Italian, German and Turkish versions

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