

Dental Wings team releases version 3.5.0.26256 of DWOS. It is an upgrade version for the iSeries, 3Series, 5Series and 7Series as well as standalone software. The upgrade brings the following new features and improvements:

### For clinics

1. iSeries -digital impression- Wizard

### For scanning

2. Default setting for hi-resolution arch scan (User preferences)

### For prosthesis design

3. Crown & Bridge Wizard
4. Introduction of the Full Virtual wax-up
5. Hollow pontic
6. Improved automatic proposition
7. Retention beads new shapes
8. New attachment kits from Bredent - free download\*
9. Introduction of attachments for bar designs from Rhein83 - free download \*
10. Introduction of "free attachments" (to add cylinder-shaped supports) and "negative attachments" (to create holes) - free download\*
11. New integrated anatomy libraries (DW Ultra and Vita's Vitapan and VPH) and more to download\* from Candulor
12. Maximum body limitation for pre-milled custom abutment design
13. Turn on/off "Automatic minimum thickness correction" in real time, when using the waxing tool on overlays
14. New insertion axis window
15. Angulated abutment screw
16. Articulator AT compatible

## For partial design

- 17. Angulated partial
- 18. Improved merging algorithms

## Model Builder

- 19. Improved stability of Model Builder
- 20. New model types: Hollow and Offset
- 21. Choice of file formats to export (Wieland, Envisiontec, Milling) according to selected material

## For guided surgery

- 22. Complete integration of coDiagnostiX
- 23. Import of specific *xorders* from coDiagnostiX

## Networking & general features

- 24. MMS is fully integrated into DWOS
- 25. Mayka Dental output available
- 26. Upgraded version (2.5) of the output file format *xorder*
- 27. Inbox local connection
- 28. Naming options for exported orders (User preferences)

\* All downloads at : <http://download.dwos.com>