

Release Notes v. 2023.1.1*

What we've fixed

We fixed the bugs below to make your CAD experience better.

Slowness in client station

Potential slowness on distributed systems have been resolved.

Issue merging crown on implant

An issue when merging a screw-retained Crown in combination with a Straumann material has been resolved. The resulting screw channel thickness is calculated according to the material specifications and can be manufactured correctly.

Gingiva scan missing

When a Synergy case is scanned with a Medit scanner the gingiva scan now appears in coDiagnostiX.

Misaligned abutments

Manual adjustment is no longer required before milling Createch-based abutment when using STL manufacturing output.

Straumann Centralized Manufacturing

Issues that could have caused delays or project cancellation in Straumann Centralized Production have been resolved.

Online Help

Issues accessing online help from the client application have been resolved.

Order creation

The issue preventing users from creating new orders when using the software in French and Spanish has been resolved.

Order Management

The issue that prevented the system from exporting manufacturing files for specific indications (Partial, Orthodontic Archiving, and Wax-up) across all language settings has now been resolved.

*This version is compatible with Dental Wings standalone software licenses and the following lab scanners:
7Series, 3Series, and Medit T500, T310, T510, and T710.